



TG16-WS Benchtop High Speed Centrifuge

TG16-WS benchtop high speed centrifuge is a classic model launched by Xiangyi/Cence for many years. Integrated with a digital display for indication of speed, time and RCF. Mature technology, stable operation and its excellent performance is deeply trusted by users. It is widely used in biology, materials science and other field.

Product advantages

Order NO	0203158000
Max Speed	16000r/min
Max RCF	21532xg
Max Capacity	6×50ml
Speed Accuracy	±30 r/min
Time Setting Range	1min ~ 99min
Noise	≤ 62dB (A)
Power Supply	AC 220±22V 50/60Hz
Power Consumption	400W
Dimension (WxDxH)	330×390×300(mm)
Packing Dimension (WxDxH)	395×650×340(mm)
Net Weight	25kg

Rotor parameters

<p>No:1</p>  <p>Rotor Type: Angle rotor Max Capacity: 12×1.5ml Max Speed: 16000r/min Max RCF: 17800xg Order NO: 0301805030 Adapter: 0.2ml, 0.5ml</p>	<p>No:2</p>  <p>Rotor Type: Angle rotor Max Capacity: 12×5ml Max Speed: 13000r/min Max RCF: 11903xg Order NO: 0301891100 Adapter: 0.5ml, 1.5ml</p>	<p>No:3</p>  <p>Rotor Type: Angle rotor Max Capacity: 12×10ml Max Speed: 12000r/min Max RCF: 14800xg Order NO: 0301805300 Adapter: 1.5ml, 5ml</p>
<p>No:4</p>  <p>Rotor Type: Angle rotor Max Capacity: 24×1.5/2.2ml Max Speed: 15000r/min Max RCF: 21532xg Order NO: 0301805510 Adapter: 0.2ml, 0.5ml</p>	<p>No:6</p>  <p>Rotor Type: Angle rotor Max Capacity: 6×30ml Max Speed: 13500r/min Max RCF: 18720xg Order NO: 0301805200 Adapter: 5ml, 10ml</p>	<p>No:7</p>  <p>Rotor Type: Angle rotor Max Capacity: 6×50ml Max Speed: 11000r/min Max RCF: 12840xg Order NO: 0301805400 Adapter: 1.5ml, 5ml, 10ml, 15ml, 20ml, 25ml, 30ml</p>