



## TDZ4K Benchtop Low Speed Centrifuge

TDZ4K benchtop low speed centrifuge is equipped with four angle rotors: 18x10ml, 12x20ml, 4x50ml 24x10ml and one swing rotor 8x10ml . It is widely used in most of the hospitals and universities by its high performance.

### Product advantages

- ◆Microprocessor control, DC brush-less motor.
- ◆Touch panel and LCD display.
- ◆Special damping structure to reduce vibration.
- ◆Integrated molding of housing.
- ◆Electronic door lock, centrifuge cannot work if the door is open and the door cannot be opened when it is working.
- ◆Countdown time displayed by seconds if the centrifugation time is less than 1min.

### Technical parameter

Order NO	0205228000
Max Speed	4000r/min
Max RCF	2200xg
Max Capacity	12×20ml
Speed Accuracy	±30r/min
Time Setting Range	1min~99min
Noise	≤ 62dB (A)
Power Supply	AC 220±22V 50/60Hz
Power Consumption	100W
Dimension (WxDxH)	456×330×266(mm)
Packing Dimension (WxDxH)	568×423×300 (mm)
Net Weight	12kg

### Rotor parameters



Rotor Type: Angle rotor  
Max Speed: 4000r/min  
Max RCF: 2200xg  
Max Capacity: 12×20ml  
Order NO: 0301808500  
Adapter: 5ml, 7ml Blood collection tube



Rotor Type: Angle rotor  
Max Speed: 4000r/min  
Max RCF: 2200xg  
Max Capacity: 18×10ml  
Order NO: 0301808400  
Adapter: 5ml, 7ml Blood collection tube



Rotor Type: Angle rotor  
Max Speed: 4000r/min  
Max RCF: 2200xg  
Max Capacity: 24×10ml  
Order NO: 0301808300  
Adapter: 5ml, 7ml Blood collection tube



Rotor Type: Angle rotor  
Max Speed: 4000r/min  
Max RCF: 2200xg  
Max Capacity: 4×50ml  
Order NO: 0301808600