



TDZ5-WS Benchtop Low Speed Centrifuge

Centrifuge TDZ5-WS is a hot-selling classic product sold by Xiangyi/Cence for many years. The sturdy metal casing is safe and reliable. Great performance, high quality and other facts make this centrifuge popular more than 20 years. Equipped with a variety of brackets and adapters, it is an ideal choice for clinical testing and laboratory analysis.



Product advantages


















- ◆ Various brackets and adapters are available.
- ◆ Microprocessor control, DC brush-less motor.
- ◆ RCF value can be calculated automatically.
- ◆ Special damping structure to reduce vibration.
- ◆ Countdown time displayed by seconds if the centrifugation time is less than 1min.








Technical parameter

Order NO	0205174000
Max Speed	5000r/min
Max RCF	4390xg
Max Capacity	4x250ml
Speed Accuracy	±30r/min
Time Setting Range	1min ~ 99min
Noise	≤ 65dB (A)
Power Supply	AC 220±22V 50/60Hz
Power Consumption	400W
Dimension (WxDxH)	450×545×380 (mm)
Packing Dimension (WxDxH)	547×727×432 (mm)
Net Weight	40kg

 Rotor parameters

<p>No:3</p>  <p>Max Speed: 4000r/min Max RCF: 2200xg Max Capacity: 2x2x96well Order NO: 0302853800</p>	 <p>Microplate Order NO: 0303853810</p>	<p>2x6 ELISA microplate</p>	2x6x96	Can be customized
		<p>2x2 Deep-well plate</p>	2x2x96	
		<p>2x1 Filter plate reagent kit</p>	2x1x96	

Rotor	Bucket	Adaptor	Capacity/Order NO	Adaptor	Capacity/Order NO
<p>No:4</p>  <p>Rotor Type: Swing Rotor Max Capacity: 4x250ml Max Speed: 4000r/min Max RCF: 3000xg Order NO: 0302853600</p>	 <p>Round Bucket Order NO: 0303891028</p>		4x1x100ml 0304891098		4x1x50ml (conical/round) 0304891109 (conical) 0304891095 (round)
			4x4x20ml 0304891101		4x4x15ml (conical) 0304891108
			4x6x15ml 0304891106		4x7x10ml 0304853604
			4x10x5ml 0304853605		4x8x1.5ml 0301891124
		Can be customized			
	 <p>Rectangular Bucket Order NO: 0303853615</p>		4x1x250ml 0304891113		4x4x50ml (conical) 0304891121
			4x7x15ml (conical) 0304891117		4x9x10ml 0304891122
			4x9x5ml 0304891196		4x15x1.5ml 0304853622

<p>No:1</p>  <p>Max Capacity: 4x100ml Max Speed: 5000r/min Max RCF: 4390xg Order NO: 0307853591</p>	 <p>Brackets Order NO: 0303853580 Capacity: 4x100ml</p>	<p>Adapter Capacity: 4x50ml Order NO: 0304891064</p>	<p>Adapter Capacity: 4x15ml (conical) Order NO: 0304891127</p>
		<p>Adapter Capacity: 4x50ml (conical) Order NO: 0304891126</p>	<p>Adapter Capacity: 4x3x5ml Vacuum tube Order NO: 0304853510</p>
		<p>Adapter Capacity: 4x30ml Order NO: 0304891125</p>	<p>Adapter Capacity: 4x2x10ml Order NO: 0304853511</p>
		<p>Adapter Capacity: 4x20ml Order NO: 0304853507</p>	<p>Adapter Capacity: 4x4 Accretion tube Order NO: 0304853512</p>
		<p>Adapter Capacity: 4x15ml Order NO: 0304891128</p>	<p>Adapter Capacity: 4x3x1.5ml Order NO: 0304853533</p>
<p>No:2</p>  <p>Brackets Order NO: 0303853545 Capacity: 4x2x100ml</p>  <p>Max Speed: 4000r/min Max RCF: 2810xg Order NO: 0307853591</p>	 <p>Brackets Order NO: 0303853545 Capacity: 4x2x100ml</p>	<p>Adapter Capacity: 4x2x50ml Order NO: 0304891064</p>	<p>Adapter Capacity: 4x2x15ml (conical) Order NO: 0304891127</p>
		<p>Adapter Capacity: 4x2x50ml (conical) Order NO: 0304891126</p>	<p>Adapter Capacity: 4x6x5ml Vacuum tube Order NO: 0304853510</p>
		<p>Adapter Capacity: 4x2x30ml Order NO: 0304891125</p>	<p>Adapter Capacity: 4x4x10ml Order NO: 0304853511</p>
		<p>Adapter Capacity: 4x2x20ml Order NO: 0304853507</p>	<p>Adapter Capacity: 8x4 Accretion tube Order NO: 0304853512</p>
	<p>Adapter Capacity: 4x2x15ml Order NO: 0304891128</p>	<p>Adapter Capacity: 8x3x1.5ml Order NO: 0304853533</p>	
	 <p>Brackets Capacity: 4x8x15ml Order NO: 0303853560</p>	<p>Adapter Capacity: 10ml Order NO: 0304853530</p>	<p>Adapter Capacity: 5ml Order NO: 0304853532</p>
		<p>Adapter Capacity: 7ml Order NO: 0304853531</p>	<p>Can be customized</p>
 <p>Brackets Capacity: 4x2x50ml Order NO: 0303853540</p>	<p>Can be customized</p>		