



CH16R Benchtop Low Speed Refrigerated Centrifuge (Specially for Blood Collection Vehicle)

CH16R refrigerated centrifuge designed for blood collection vehicle. Swing rotor 16x15ml and adapters, accommodate to vacuum blood collection tube, nucleic acid centrifuge tube and 10ml 15ml tubes. It's specially suitable for nucleic acid centrifugation detection. Stainless steel Chamber, multiple protection, safe and reliable. The compressor is fully enclosed type and rotary type. CH16R is really a good choice for blood collection vehicle.

Product advantages

- ◆ With unique sucker type foot, it can be sucked on the table tightly.
- ◆ With 5 inch HD touch screen control, easy- operate and intuitive. The setting parameters and running parameters at a glance.
- ◆ Multiple adapters available, it's a compact, versatile and indispensable general purpose centrifuge.
- ◆ Environmental friendly compressor is adopted, which has good refrigeration effect and no pollution.
- ◆ Multiple protection, include over-speed, over-temperature, lock and so on, to ensure user's safety.

Technical parameter

Order NO	0204260000
Max Speed	5000r/min
Max RCF	3555xg
Max Capacity	4x4x15ml or nucleic acid tube
Speed Accuracy	±30 r/min
Time Setting Range	1min~99min
Temperature Setting Range	-20°C ~ 40°C
Temperature Accuracy	±1°C
Compressor Unit	High performance compressor unit , Environment-friendly refrigerant (134a)
Noise	≤ 65dB (A)
Power Supply	AC 220±22V 50Hz
Power Consumption	500W
Dimension (WxDxH)	415×625×320(mm)
Packing Dimension (WxDxH)	960×530×520(mm)
Net Weight	55kg

Rotor parameters

No:1



Rotor Type: Angle rotor
Max Speed: 5000r/min
Max RCF: 3555xg
Max Capacity: 12x15ml
Order NO: 0301808740

No:2



Rotor Type: Swing Rotor
Max Speed: 4200r/min
Max RCF: 2958xg
Max Capacity: 4x4x15ml
Order NO: 0302855740
(Brackets not included)



Brackets
Capacity: 4x50ml
Order NO: 0303855210
Adapter: 15ml (conical)
Order NO: 0304891103



Brackets
Capacity: 4x4x15ml
Order NO: 0303855745
Adapter: 1.5/5/10ml
0304854806 1.5ml
0304854808 5ml
0304855729 10ml



Brackets
Capacity: 4x6x5ml
Order NO: 0303855230
Adapter: 1.5/0.5ml
0304854807 1.5ml
0303854813 0.5ml