



## WTL-4K/6K/10K

### Micro Low Speed/ High Speed Centrifuge

WTL-4K/6K/10K is a small but powerful mini centrifuge for simple and quick liquid separation. With an 6 place 1.5/2mL tube rotor as well as a 2 place PCR 8-strip tube rotor, WTL mini centrifuge can accommodate 0.2mL 8 strip tubes or 16 individual PCR tubes. With the included adapter, 0.2mL and 0.5mL tubes can also be accommodated in the 6 place tube rotor. This centrifuge provides quiet, safe, and easy operation for your everyday quick spin-down applications.

### Product advantages

- ◆ Voiceless technology, running with low noise
- ◆ Equipped with brushless motor, service life is long.
- ◆ Durable and safe for users.
- ◆ Equipped with an 6 place 1.5/2mL tube rotor as well as a 2 place PCR 8-strip tube rotor.
- ◆ AC100 to 242V, global application.
- ◆ WTL mini centrifuge takes simple one-step operation to a new level. Simply place your sample in the centrifuge, close the lid, and it will begin to rotate without the press of a button.



### Technical parameter

| Model                     | WTL-4K                         | WTL-6K                         | WTL-10K                        |
|---------------------------|--------------------------------|--------------------------------|--------------------------------|
| Order NO                  | 0205220100                     | 0205220200                     | 0203220300                     |
| Max Speed                 | 4000r/min                      | 6000r/min                      | 10000r/min                     |
| Max RCF                   | 900xg                          | 2000xg                         | 5000xg                         |
| Max Capacity              | 6x1.5ml/0.5ml/0.2ml, 2x8x0.2ml | 6x1.5ml/0.5ml/0.2ml, 2x8x0.2ml | 6x1.5ml/0.5ml/0.2ml, 2x8x0.2ml |
| Noise                     | ≤53dB                          | ≤53dB                          | ≤53dB                          |
| Power Supply              | AC100~242V                     | AC100~242V                     | AC100~242V                     |
| Power Consumption         | 30W                            | 30W                            | 30W                            |
| Dimension (WxDxH)         | 160x200x123(mm)                | 160x200x123(mm)                | 160x200x123(mm)                |
| Packing Dimension (WxDxH) | 205x265x160(mm)                | 205x265x160(mm)                | 205x265x160(mm)                |
| Net Weight                | 1.2kg                          | 1.2kg                          | 1.2kg                          |